


Wellman Plastics Recycling Site – Johnsonville, SC

(2 million square feet under roof)



Our PCR Nylon Raw Material Flow



RM to FG Viscosity (dl/g) Comparison

Material	Average	Standard Deviation
Virgin Nylon 66	1.54	0.09
EcoLon® 2100-BK1 (100% PCR) – RM	1.50	0.03
EcoLon® 2100-BK1 (100% PCR) – FG	1.35	0.05



EcoLon® consistency is better than Virgin Nylon



WELLMAN
Engineering Resins